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

Electrode Usage Guide


1. To avoid dermatological infection and for hygiene purpose, please use dedicated pad for each person.
2. Please keep the water-based gel on pad clean and avoid direct touching.
3. Place the pad on clear plastic cover to avoid contamination and extend usage life.
4. Please keep pads under shade and avoid high humidity and direct sunlight. Recommended storage condition: Temperature <28°C and humidity 50±5%
5. The average pad usage life is around 30~40 times (roughly one to two weeks) depending on personal hygiene, skin condition and how many times use per day.

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⚠ Other Considerations


1. Results from UGYM Pro may differ from person to person.
2. Please follow doctor's instructions when using UGYM Pro if you are currently undergoing physical therapy. Visit a doctor from time to time to better understand the recovery progress.
3. If any skin allergy or irritation is caused by the pad gel, stop using UGYM Pro immediately and consult a doctor.

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▼ Recommended pad placement

The actual feeling of stimulation intensity relies on pad locations, physical condition, size of electrodes, distance between each electrode and personal sensitivity.

Please drink 500cc of warm water after each massage session to keep muscle cells hydrated during recovery and speed up metabolism of waste product.

Quick Guide

1. Thoroughly wash and dry the skin, snap the black or red lead to the electrode clip button, slowly pull electrode off the clear plastic sheet by grabbing its ear, apply electrodes flat to the skin. Make sure no electrodes are overlapping. Uneven electrodes will affect the desired stimulation.
2. UGYM Pro has 4 massage modes to choose from. Select the mode that fits you best and raise the intensity level where you find most comfortable. Be careful not to raise the intensity level too fast as overloading intensity may cause muscle cramps or tear.
3. Every massage session lasts for 30 minutes and will stop automatically. If you change the massage mode during previous session, the timer will restart for 30 minutes.
4. If you want to remove electrodes during massage session, remember to press reset button or switch to the lowest intensity. (Reset button is located next to “+” button)
5. Keep in mind you can use your smartphone to control the device remotely. (Bluetooth version only)

Troubleshooting

Device has suddenly stopped during session.

1. Check to see if the battery LED display is blinking in orange. The device may need to be charged.
2. Check to see if electrodes are falling off the skin; one of the leads has got disconnected; or if the lead cord jack got disconnected from the device.
3. Check the lead cord for cracks, tears or other possible damages.
4. Every massage mode lasts for 30 minutes and the device will automatically turn off after 60 seconds with no further action.
5. If your device is overheating or malfunctioning, it will automatically power off.

How to reset the device if it's malfunctioning?

1. Leave it for a day and it will become normal again.
2. Let the battery drain and recharge the device again. (May take couple days)
3. Charge the device for 10 minutes and unplug the micro USB cable.

⚠ Embedded 3200mA Lithium Battery

UGYM has a battery life up to 8 hours. Depending on the program and intensity selected, the battery life may vary from person to person. Through the USB port, UGYM can also be used as an external battery to charge other electronics. When the device is fully charged the LED will display green. As the battery life diminishes, it will blink orange. All functions will shut off during charging.

*** Attention

- Do not use U-GYM in temperatures above 37°C.
- Do not put U-GYM in or around fire.

About Electrodes (Pads)

1. To avoid dermatological infection and for hygiene purpose, please use dedicated pad for each person.
2. Please keep the water-based gel on pad clean and avoid direct touching.
3. Place the pad on clear plastic cover to avoid contamination and extend life after use.
4. Please keep pads under shade and avoid high humidity and direct sunlight. Recommended storage condition: Temperature <28°C and humidity 50±5%
5. The average pad usage life is around 30~40 times (roughly 3 to 4 weeks) depending on personal hygiene, skin condition and how many times use per day.

Electrodes and lead cords are consumables with limited life expectancy and no warranty provided.

UGYM Technology Corp.

No.2-1, Ln. 18, Yanshou St., Songshan Dist., Taipei City 105, Taiwan
Tel: 886-2-27476176 Email: info@myugym.com

[f](#) [i](#) [t](#) @RelaxDeepMuscle

RelaxDeepMuscle.com

U>GYM

Empowering Muscle Strength

UGYM Pro Smart Deep Tissue Muscle Massager

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


NCC Notice
According to NCC regulation for low frequency and radiation electronics section 12: Company, business and end user can not alter frequency, enhance power or change original certified design and functionality without official authorization; section 14: Low power radiation electronics cannot be used to affect flight safety or interfere with approved communication and should stop using the electronics immediately if situation appears so.

This device complies with part 15 of the FCC rules. Operation is subject to the follow Two Conditions.
(1) This device may not cause harmful interference and
(2) This device must accept any interference received including interference that may cause undesired operation.

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




How to download app
UGYM Supports iOS & Android
ugym


Just use keyword "UGYM" in app store and you should find UGYM logo. Please download and install the app in your smartphone.

2 Activez l'application

Once you open the app you will see the login screen (right). You can either use Facebook login or Airsig login.


When using Facebook to login, make sure the Facebook account is already login, the app will take you to Facebook login page and select "Use Facebook to Login" and click "Confirm" to continue. The app now will take you to the main page. (Note that if you are inactive for 5 minutes, you will be logged out.)


When using Airsig to login the first time, the prompt will lead you to setup a signature. Please think of a name at least 6 characters such as "Franky", picture that word in your head and keep your thumb on the screen and hold the phone like a wand to write the word as if you are writing on the wall with the tip edge of your phone. Once it's setup, it will take you to the main page. (Note that if you are inactive for 5 minutes, you will be logged out.)

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You will need to pair the device to update the firmware or control the device. Press the upper left corner and click on "Settings" and click on UGYM Paring. Now you see a device turned on named "U-Gym 075DBC", click on the name and you can set a nickname for it. Now the device is paired successfully. You can click on "System Update" to see if there's any available firmware to update your device.

Now click on "Mode" to go back to the main screen.


Click on one of the M1 - M4 modes and the screen will change to below. Press "+" to go up higher intensity and "-" to go down. You can see the remaining time left on the lower bar.

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
Press the upper right corner and select "Download Available Modes" and you will see UGYM005 to UGYM011 show up in the Available Mode section.


If you select one of the available modes download from our server, you can replace M1 to M4 mode on the device. If you don't like the change, you can restore the default M1-M4 back on device as well.

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Before using the Self-define mode, please make sure you click on upper corner to select "Sync Available Cycles" and you will see different cycles pops up.

Click on any of the cycles and they will show up in the row of sequence. Click on "Experience" and intensity to try your design. (Make sure your pads are already on your body and cord is plugged in to the device)

If the design is complete, press "Save" and give the mode a name. Now you can see there's additional mode on the bottom of main page.

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Click on one of the Self-design mode and you can delete or share the mode on Facebook.