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
How to check battery life?

Long press 0.3 second to turn on the device. While green or blue LED is blinking, press power button once and white LED will start to blink. Blink once indicates battery life between 10 to 39%; blink twice for 40 to 59%; blink 3 times for 60 to 79% and blink 4 times for more than 80%.

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
Operation Notice:

1. Please clean the skin using water and keep it dry; lift up pad (electrode) from rounded ear side off the clear plastic cover and place the pad to the treated area. Keep the pad flat on skin and avoid any wrinkle for better effect.
2. There are two stimulation modes to choose from. Select the desired mode and press + or - to change the intensity level. Slowly increase the intensity level and keep intensity bearable and comfortable. (Double press ⏻ to switch mode)
3. The stimulation mode is set to end automatically in 30 minutes. The timer and intensity level will reset if the mode has been switched.
4. UGYM mini will shut off automatically in 60 seconds if no further action is followed.
5. Please do not remove the pad from skin until the intensity level is switched to the lowest level or device is turned off.
6. Stimulation cannot be activated for battery lower than 10% and device cannot turn on if battery is lower than 5%.

U-GYM simulates deep massage technique from real therapist to give you full relaxation!

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

Supported Electropad


Currently we have three sizes to choose from


Size 5X5cm, 4 pads a pack

Size 5X9cm, 4 pads a pack

Size 7X13cm, 2 pads a pack

※Recommend to use 4 pads with 1 to 4 cord;
and use 2 pads with 1 to 2 cord

Notice during Massage:

The stimulation intensity is selected based on different treatment areas, size of pads, personal muscle sensitivity and individual body condition. Make sure the intensity is at your most comfortable zone.

- 1.1 hour (two cycles) of treatment for each area per day is recommended.
- 2.Please drink 500cc of warm water after each session to keep muscle cells hydrated during recovery and speed up metabolizing extra cell waste.
3. The actual feel of stimulation intensity relies on pad locations, personal physique and individual sensitivity.
4. Please stop using UGYM mini immediately if you are feeling discomfort or ill in any way during the session and consult your doctor.

Instruction for Electrode:

- 1.Avoid placing pads at heart, neck, eye, oral, genital, arteries and sensitive areas.
- 2.Avoid putting pads at rash, zit and wounded areas.
- 3.For users with skin disease history, please be careful not to use over time.
- 4.If the pads are contaminated, dried out, worn off or other factors causing weak stimulation performance, we recommend to replace with new pads. (For best stimulation result, we recommend to replace new pads every month or after 30~40 uses.)

Best Practice for Electrodes:

- 1.To avoid dermatological infection and for hygiene purpose, please use dedicated pad for each individual.
- 2.Please keep silicone gel on pad clean and avoid touching directly.
- 3.Place the pad on clear plastic cover after every use to avoid contamination and extend usage life .
- 4.Please keep pads under shade and avoid high humidity and direct sun light. (Recommended storage condition: Temperature <28 degree and humidity 50±5%.)
- 5.The average pad usage life cycle is around 30~40 times (roughly within one month) depending on use frequency, personal hygiene and skin condition.

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
Please contact us for assistance:
E-mail: info@myugym.com
Phone: 886-2-2747-6176
Hours: Monday to Friday from 10AM to 6PM

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Notice for Charging:

If orange LED light is off while charging, that indicates the device is fully charged.

Embedded 2600mAh Lithium Battery:

1. UGYM mini can be used up to 5 hours after fully charge (Usage time varies and is depending on different intensity level and sequence used)
2. The orange LED will flash slowly to indicate low battery level (Less 10%).
3. Stimulation will be disabled during battery charging. Both charging and stimulation cannot perform at the same time.


Warning:

Please do not use UGYM mini in environment over 40 degree Celsius.
Throwing the device and battery in fire is strictly prohibited.

Troubleshooting - Device stops during stimulation:

1. Please check if the LED is flashing in orange showing low battery. Please charge the device immediately if that's the case.
2. Whether electrode cord is firmly plugged in the device and all snap-on clips are firmly attached to the pads. Also, check to see if any of the electrodes are fallen off from skin.
3. Please check if electrode cord is damaged or cracked.
4. If the session has been run for 30 minutes, the device will shut off automatically.
5. If the device is overheating or acting abnormally, it will shut down automatically for safety.

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U-GYM Technology Corp.